

PRO25

Audiometric mini booth

puma
soundproofing



The smallest, lightweight and portable soundproof booth ready to be used in a few minutes.

Designed for hearing screening in audiology centres.

Available in DIY kit with reversible door and magnetic closing.

Folding table and connecting panel



All our cabins offer the best of Italian design



Product conforms with directive EEC/93/42 regarding class I medical devices

Noise reduction

PRO25 is tested by means of acoustic insulation tests at the Technical-Physics Laboratory C.S.I. S.p.A. (Milano-Italy) (ref. Standards UNI EN ISO 11957 - UNI EN ISO 717-1).

Acoustic insulation

Hertz	125	250	500	1000	2000	3150	4000	6300	8000
dB	18,0	16,0	26,7	30,8	28,8	30,0	30,9	38,0	39,0

Applications

PRO25 is recommended for hearing screening in audiology centres for those who have limited space available and need a very small soundproof booth.

PRO25 can be also used in temporary situations for occupational medicine, in companies and schools.

Configuration

Available in a single standard size.

External size: 97x97x192 cm (l x w x h)

Internal size: 83x83x180 cm (l x w x h)

Weight 175 kg.

Exterior

The exterior finish is made of a white, highly-resistant, fireproof, washable panel.

Interior

The interior finish is made of acoustic pyramid foam in anthracite grey, designed to reduce sound waves for a perfect acoustics.

Door

The door is reversible with magnetic closing. The opening (right or left) can be decided upon assembling the booth and changed in 5 minutes. Clear opening 55x168 cm (w x h).

Window

The window consists of a special double safety glass. Clear view 56x56 cm (w x h).

Floor

The floor is layered with different materials and it stands on anti-vibration adjustable feet.

Finish: technical antistatic washable carpet.

Supports for instruments

External folding table for instruments (79x50 cm) and internal hangers for earphones are included.

Connecting instruments

Cables are laid through a new universal soundproof panel without jack connection. All kinds of cables can be used without the need to be cut. Adaptable to all types of instruments.

Electrical system

Lighting is ensured by a fluorescent lamp with low energy consumption.

Power supply: 220/240 V. - 50/60 Hz.

Cable for power supply included.

Delivery

Delivered in packaging for shipping in a DIY kit.

Standards

MDD 93/42/CEE - 2007/47/CE - CEI 64/11

CEI 64/8 - 2006/95/CE.

Installation

Easy and quick to install; no specialised staff required. Just follow the instructions.

Thanks to reduced weight and small size PRO25 can be easily set up in 5 minutes and transported in a Station Wagon car.

Notes

The materials are fireproof and self-extinguishing at the origin.

Specifications are subject to change without notice.

Warranty: limited to one year from delivery date.

