



enterMed

PATIENT CHAIRS

Entermed Ergonomic Examination/Treatment Chair

- Electrical height adjustment, foot controlled
Min. seat height: 51.0 cm
Max. seat height: 81.0 cm
- 360° rotatable
- Backrest adjustable from vertical to horizontal position by gaspring system
- Headrest and armrests adjustable
- Optional hinged footrest
- Optional automatic return facility
- Standard colours black 201, pigeon blue 134, pearl grey 129

Product codes

- | | |
|----------|--------------------------------|
| 760103 | Patient chair without footrest |
| 760105 | Patient chair with footrest |
| 760105.5 | Automatic return facility |



760105



Product code:
750109

ENT Examination/Treatment Chair

Product code: 750109

- Synchronised adjustment of backrest and footrest to horizontal position
- Stylish design
- Electrical height adjustment, foot controlled
Seat height min.: 48,0 cm, max: 88,0 cm
- 180° rotatable
- Very stable
- Removable Armrest
- Easy to clean (seamless upholstery)
- Standard colours black 201, pigeon bleu 134 and pearl grey 125

Budget-friendly ENT Examination Chair

Product code: 750108

- Low cost
- Electrical height adjustment, foot controlled
Seat height min.: 53,5 cm, max: 69,5 cm
- 360° rotatable
- Backrest adjustable to 160° by means of lever
- Armrest can be turned backwards or removed
- Easy to clean (seamless upholstery)
- Standard colours black 201, pigeon bleu 134 and pearl grey 125

Product code:
750108



enterméd

DOCTOR'S STOOLS

UNIVERSUM bucket stool

- Gasspring height adjustable by means of lever
Min. seat height: 41.0 cm
Max. seat height: 54.0 cm
- 360° rotatable
- Seat inclination device
- Anatomic design
- Standard colours black 201, pigeon blue 134 and pearl grey 129

Product code
753000 UNIVERSUM bucket stool

Product code:
753000



REGULAR doctor's stool

- Gasspring height adjustable by means of lever
Min. seat height: 46.0 cm
Max. seat height: 58.0 cm
- 360° rotatable
- Adjustable backrest (height and depth)
- Seat Ø 36 cm
- Cross-type base chrome-plated
- Braked safety castors
- Standard colours black 201, pigeon blue 134 and pearl grey 129

Product code
30142 ES Regular doctor's stool



Product code:
30142 ES

Svensk distributör:
Entomed MedTech AB
Storgatan 7
233 31 Svedala
040-481140

ENTERMED B.V.
P.O. Box 236
3440 AE Woerden
The Netherlands
Tel : + 31 348 415 184
Fax : + 31 348 422 064
Internet: www.enterméd.com

Showroom address:
Laan van Overvliet 7
3461 HE Linschoten

